

MEDIVENT EXCEL

Electronic Anesthesia Ventilator

MEDISYS[®]

MEDIVENT EXCEL

Electronic Anesthesia Ventilator

Product Catalogue



MEDIVENT EXCEL
Electronic Anesthesia Ventilator



SPECIFICATIONS

- CLASSIFICATION
- MODES
- TIDAL VOLUME
- FREQUENCY
- INSPIRATORY TIME
- I/E RATIO
- INSPIRATORY HOLD
- APNEA BACKUP
- ALARMS
- SAFETY VALVE
- DRIVING SOURCE
- POWER CONSUMPTION
- WEIGHT
- DIMENSIONS

Flow Generator Type

Cmv, Cmv+assist & Simv

50 To 1200ml

6 To 40bpm (1 To 40 Bpm In Simv)

0.4 To 3 Seconds

1:1 To 1:5

0 To 2 Seconds

Every 4 Secs In Case Of Apnea.

Patient Disconnection, High/low Pressure

10 To 60 Cm H₂O

220v Ac

60 W

12 Kgs.

410x5800x750mm

OUR PRESENCE

