

SOFT TISSUE DIODE LASER

The Premier & Most Complete Diode Laser



Power & Versatility in a Small, Portable Package

940 nm is a superior diode wavelength

- ④ Faster treatment and healing
- ④ Clinical versatility for complete patient care
- ④ Portable for easy transport
- ④ ComfortPulse™ modes for the ultimate patient comfort

Intelligent Procedure - Designed to Maximize Efficiency

- ④ Easy touch-screen interface and pre-sets to shorten your learning curve
- ④ Most common and productive procedures already programmed
- ④ Totally wireless laser for easy portability between operatories

EPIC X Technical Specifications	
Laser Classification	IV
Medium	InGaAsP Semi-conductor diode
Wavelength	940 ± 10nm
Max Output Power	10 W
Pulse Repetition Rate	upto 20 KHz
Pulse Duration Rate	0.01 ms - 20 ms
Pulse Interval	0.04 ms - 20 ms
Power Modes	Continuous, Pulse Modulation
Aiming Beam	Laser Diode, max 1mW, 635nm
Fiber Tips Diameter	200 µm, 300 µm, 400 µm
Dimensions	5.7in (W) x 4.47in (H) x 6.5in (L) (14.5cm x 11.2cm x 16.5cm)
Weight	2.5lb (1kg)



FDA CE ISO 13485:2012 + AC:2012

Surgery, Pain Relief & Whitening At Your Fingertips




Create Bright Smiles

EPIC X offers the fastest in-office whitening procedure available, with dramatic shade change in 20 minutes of chair time.

Provide Pain Relief

Treat TMJ & other myofacial disorder with patented handpiece attachment to simplify the procedures.

Maximum Performance, Minimum Chair Time



Ver 12/2022