

Vital Signs Monitor

Flight[®]
ONE

FEATURES

- Compact and flexible appearance.
- Simple and convenient operation, which is applicable for community medical treatment and family etc.
- 4.3 Inch color TFT display and highlight LED display adjustable and audible and visual alarms.
- Built-in rechargeable lithium battery for uninterrupted
- Monitoring
- NIBP, SpO2, Temperature

NIBP

- Method: Oscillometry
- Measure mode: Manual, Auto, STAT
- Unit: mmHg, kPa
- Measure and alarm range:
- Adult Mode
 - SYS 40~280 mmHg
 - DIA 10~225 mmHg
 - MEAN 20~240 mmHg
- Pediatric Mode
 - SYS 40~220 mmHg
 - DIA 10~160 mmHg
 - MEAN 20~170 mmHg
- Neonatal Mode
 - SYS 40~135 mmHg
 - DIA 10~100 mmHg
 - MEAN 20~110 mmHg
 - Resolution: 1mmHg
 - Accuracy: ±5mmHg

Standard Parameters:

- NIBP, SPO2, PR, Temp.



Technique Specification

- SPO2
- Measuring Range: 0~100 %
- Resolution: 1%
- Accuracy: 70%~100% (±2%)
- Pulse Rate: 20-250 BPM
- Resolution: 1 BPM
- Accuracy: ±3 BPM
- Temperature probe

Standard Accessories:

No.	Item	Qty
1	Main Unit	1
2	Adult SpO2 probe	1
3	Adult NIBP cuff	1
4	NIBP extension tube	1
5	Temperature probe	1
6	Power Cable	1
7	User Manual	1

Packing:

- Single package size: 35* 25* 26cm
- Single gross weight 3KG
- 8 Unit per carton, package size: 72.5*53
- *56.5 cm, total gross weight 24KG