



RAD^X | HFX 5



High Frequency Mobile X-ray


5 kW with APR

HFX 5 mobile X-rays is based on high frequency and high voltage X-ray generator, which ensure superior image quality, with low skin dose to protect the patient effectively.




● Useful Output ● Skin doses ● X-ray not useful

HF X-ray generator use the frequency-inverter-rectifier technology, to produce near constant potential



Affordable, Digital Upgrade

DR Digital X-ray Detector *PLUG & PLAY*



- Portable digital flat panel detector
- 17" x 17" active area
- Long life battery
- High SNR and Low Noise at lower dose
- Faster image acquisition in less than 3 second
- Large on-board storage memory
- Add on reduced patient dosage
- Wire / Wireless image transfer
- Light weight with enhanced durability
- Reliable high-performance imaging
- High sensitive AED function
- Easy Grip and removable handle

Flat Panel X-ray Detector

(optional)



Dual Battery



Status Monitoring OLED



Low Dose IGZO TFT



Enhanced AED®



Extended Battery Life



Durable Design



Dust and Water Proof

Realtime Image Display

