

Portable Ventilator

AIRLOG 3

Standard Configuration

Main Unit, Ventilation Hose, Patient valve, Mask

Features

- IPPV, Assist & Manual
- TV range: 75-2000ml
- Breathing rate: 5-40 BPM
- Inspiration : Expiration 1 : 1.67
- Minute volume (MV): 3-20 LPM
- Pressure limit: 20-60 mbar
- Trigger at less than -2 mbar
- PEEP range: 5-20 mbars
- Alarm/warning once patient stop breathing
- Gross weight not over 1.3kg (just only the machine, not include tube and etc.)
- IP class: IPX4
- Battery working time: 10 hours

Application

- Wild rescue, vehicle transport, long-distance transport, hospital transport



CE