

360° ROTATE CONTRA ANGLE



360°

The contra angle can be 360° rotated without taking of to make it easy to use during treatment



135° C High temperature and high pressure sterilization

MINI HANDPIECE
WITH 16:1 ELECTRIC CONTRA ANGLE



Mini head



Normal head

 **MEDIKABAZAAR**

Healix®

Wireless Endomotor

- ✓ Mini head Hand Piece
- ✓ 9 Custom Programs
- ✓ 360° Rotary Angle
- ✓ Colorful OLED Screen



Marketed By:

 **MEDIKABAZAAR**

Boston Ivy Healthcare Solutions Private Limited

3rd Floor, Second Avenue, Atul Projects,
Subhash Nagar, Andheri East, Mumbai,
Maharashtra 400093.

 +91 794805 8625

 support@medikabazaar.com

 www.medikabazaar.com

Features:

1. Ergonomic design body, light weight and easy to carry
2. With mini push button 16: 1 contra angle(left and right reciprocating), can withstand 135°C high temperature and high pressure sterilization.
3. 10 kinds of procedures(memory set) for selection.

Technical:

- Battery voltage: 3.7V
- DC Battery capacity: 800mAh
- Speed range: 120-650rpm
- Torque range: 0.6-4.0N.cm
- Rotation angle: 30-370°

16:1 Angle Handpiece

Power key/Motor Switch

Display Screen

Item Select Key

+Key

-Key

5V in

9 Memory Modes

Up to 10 different programs can be set to adapt to different clinical operation requirements

M9

120rpm



REV 4.0N.cm



Press-to-remove root canal files
Safe and fast



Comfortable grip

no fatigue after long-term use



Reciprocation mode

3 working modes: **FWD**, **REV**, **REC**



Motor

REV

FWD - Forward rotation(Automatic reversion of torque)

REV - Reverse rotation,

REC - Reciprocating motion

(Angle is adjustable and combination freely)