

Tropicalgin

multiPro

Chromatic alginate for versatile use (Fast and Extra Fast).



PRODUCT BENEFITS

- **Chromatic:** provides the practitioner with a visual guide during material mixing, working and positioning in the oral cavity
- **Different setting times (Fast/Extra Fast)*:** to satisfy all dentists' requirements
- **Mango scent:** favours patient comfort



MANGO

The chromatic variation provides the practitioner with a visual guide during the different phases:

- red during mixing
- orange during the working phase
- yellow during positioning in the patient's mouth

Clinical times**

	Fast	Extra Fast
Working time (including mixing time)	1' 35"	1' 15"
Time in mouth	1' 00"	25"
Setting time	2' 35"	1' 40"

* Product not available on all markets

** Clinical use times are intended as from the start of mixing with deionised water at 23°C / 73°F.