



MEDIKABAZAAR



NAVAKAR
DENTAL SUPPLY CO.

INDEX

	Page No.		Page No.
1 MTA-PUTTY	3	28 LIDO-PULP	17
2 BIOSEAL - ACTIVE	3	29 ALVEOLAR	17
3 BIO - REZ	4	30 ND-TEMPFILL	18
4 DUO-CORE	5	31 APEX-SEAL	18
5 RESIFIX – DUO	5	32 F-VARNISH	19
6 INTAGLIO	6	33 N-SDF	19
7 INTAGLIO – REFILL PACK	6	34 ANTIFOG SOLUTION	20
8 PERPLEX FLOW	7	35 ROTO SPRAY	20
9 TETRA FLOW	7	36 CARIS DETECTOR	21
10 ALL IN ONE FLOW	8	37 ENDO-SOL	21
11 PIT AND FISSURE SEALANT	8	38 DECENS-PF	22
12 MICRO LINER	9	39 FORMO KRESOL	22
13 THIXOCAL-LC	9	40 UNILIK	23
14 PROSEAL	10	41 PREP-AQUA	23
15 CALSILLINER	10	42 EUGENOLL	24
16 α BOND	11	43 COPAL VARNISH	24
17 α BOND 007	11	44 THIXOCAL	25
18 MAXXWHITEOFFICE	12	45 THIXO-PREP	25
19 MAXXWHITEOFFICEPLUS	12	46 STYPTIC-GEL	26
20 MTA	13	47 IODOPEX	26
21 MTA-JET	13	48 BRIGHT-SMILE	27
22 AFFIX CEM ONE	14	49 ND-SHINE	27
23 AFFIX CEM TWO	14	50 THIXOETCH	28
24 NATURAL – CARE	15	51 THIXO ETCH JUMBO	28
25 PORCELAIN PREP KIT	15	52 RVG SLEEVES	29
26 ORAGEL	16	53 ZINC –OXIDE	29
27 SPARKLE	16		

MTA - PUTTY



packaging

2 x 1gm Syringe

Premixed Injectable MTA Putty
Contains : Calcium Silicates, Calcium aluminates, Calcium alumino-ferrite, Calcium Sulfates, Radiopacifier

INDICATIONS FOR US :

- Root canal filling.
- Furcation perforation repair.
- Resorption repair.
- Pulp capping.
- Pulpotomy.

Bioseal - Active

(Premixed Injectable bio ceramic sealer)

Contains : Calcium Silicates, Calcium aluminates, Calcium alumino-ferrite, Calcium Sulfates, radiopacifier



INDICATIONS FOR US :

- Bioseal - Active is developed for root canal filling and root perforation repair, effectively preventing periapical tissue irritation as well as secondary infection.
- Simple and fast canal filling.
- Strong antibacterial effects.
- Promotion of hard tissue formation- No postoperative pain after overfilling. - Complete apical sealing.

packaging

1 x 3 gm Syringe
& dispensing tips

Bio - Rez

Self Etch Self Adhesive Resin Cement



INDICATIONS FOR US :

- Self-adhesive resin cements are designed to simplify the placement of indirect restorations by eliminating the need for separate etchants and primers.
- Dual cured and available in different shades to match the clinical situation at hand, these cements are most often for use with ceramic or metal crowns, bridges, inlays, onlays and endodontic posts.
- Most self-adhesive cements provide good bond strengths and can save clinicians time because fewer steps are required.
- As with any cement system, be sure to use the self-adhesive resin you choose according to the manufacturer instructions for placement and curing..

MADE IN INDIA
MADE FOR WORLD

Duo – Core

Dual cure auto mix nano filled Radiopaque composite core Buildup Material



Core build-up materials are used when there is not enough healthy tooth left for proper crown placement. Building up the core of the tooth with resin or plastic composites allows for better retention of the crown and reduced sensitivity while strengthening the tooth itself.

packaging

8 gm syringe

Working Time 90 Sec.
Curing Time 4 to 5 min.

Resifix - Duo

Dual cure Resin Modified Radiopaque Glass Ionomer Luting Cement



packaging

8 gm syringe


INDICATIONS FOR US :

- Permanent cementation of crown and bridge.
- Permanent cementation of endodontic posts made of metal, reinforced composite.
- Permanent cementation of implant retained crowns.
- 0.2mm film thickness

Working Time 1.30 min.
Setting Time 4 to 5 min.

INTAGLIO COMPOSITE KIT

Light curing Nanobase Dental Restorative Material



Benefits & features

- High compressive strength, highly radiopaque.
- low shrinkage.
- Compressive Strength 395 Mpa.
- Curing time 20 sec - 40sec depending upon shade depth.
- Used for : class I, II, III, IV, V Restoration.
- Direct Aesthetic Veneers for Masking & Anomalies.
- Composite & Porcelain Repairs.

INTAGLIO – REFILL PACK



Advantages

- Radiopaque
- Highly realistic
- Brilliant shine
- Low curing shrinkage
- 7 shades available
- Excellent malleability

compressive
strength
395 Mpa.

packaging

1 syringe 4 gm

INDICATIONS FOR US :

- Direct restorations of class I to V cavities 7 shades available.
- Restoration of deciduous teeth.



Packaging

2 x 2 gm Syringe

Tetra Flow

Lightcured Nano hybride flowable composite

Contains : Methacrylate monomer resins

INDICATIONS FOR US :

- Class V restorations.
- Sealing margins.
- Repairing pitted areas in hybrid restorations.
- Small occlusal and lingual repairs.
- Sealing gingival margins in class II and III restorations.
- In deep class V restorations to block out shine through prior to overlay to tetra flow nenofit.

Perplex Flow

Self Adhesive Light cured flowable composite

Contains : Methacrylate monomer resins

INDICATIONS FOR US :

- Class V restorations.
- Sealing margins.
- Repairing pitted areas in hybrid restorations.
- Small occlusal and lingual repairs.
- Sealing gingival margins in class II and III restorations.



packaging

2 x 1.5 gm Syringe

All In One Flow

Lightcure self etch self adhesive flowable composite

Contains : Methacrylate monomer resins & barium glass fillers

INDICATIONS FOR US :

- The composite is immediately ready for use because the preparatory steps of etching, priming and bonding are contained in Constic
- t cementation of crown and bridge.
- Permanent cementation of endodontic posts made of metal, reinforced composite.
- Permanent cementation of implant retained crowns.

packaging

2 x 1.5 gm Syringe

Pit & Fissure Sealant

Methacrylate monomer resins & fillers

INDICATIONS FOR US :

- Sealing of teeth with sound but valnerable pit & fissure sites.
- Sealing of teeth with minimal pit & fissure caries.
- Acts as barrier to plaque and bacteria.
- Can also be used to seal margins of composite restorations.

packaging

2 x 1.5 gm Syringe

Thixocal - LC

Lightcured radiopaque calcium hydroxide paste

Contains : Methacrylate active monomer, calcium hydroxide, anti oxide, radiopaque material.

INDICATIONS FOR USE :

- Used as cavity liner and pulp capping.
- Metal opifier in prosthesis repair.



packaging

2 x 1.5 gm Syringe

Micro Liner

Light cure Glass Ionomer composite Liner

Contains : Methacrylate monomer resins
Glass Ionomer powder

INDICATIONS FOR USE :

- Forms a strong bond to dentin, sealing the tooth to reduce the incidence of microleakage, which helps reduce post-operative sensitivity.
- Produces zones of inhibition against bacteria commonly found in the oral cavity in-vitro.



packaging

2 x 1.5 gm Syringe



packaging

2 x 1.5 gm Syringe

Calsil Liner

Resin Modified Calcium Silicate Filled Liner

Contains : Calcium Silicate

INDICATIONS FOR USE :

- Calsil liner precise placement allows its use in all deep cavity preparations.
- The light-cured set permits immediate placement and condensation of the restorative material.
- Its proprietary formulation allows for a command set with a light curing unit while maintaining ease of placement due to thixotropic properties.

Proseal


Lightcure Protective & Glazing Surface Sealant


INDICATIONS FOR USE :

- Proseal is the only light-cured, protective clear coating formulated with adhesive monomer and nano fillers. One thin coat protects margins, prevents staining and provides a high gloss (comparable to intense polishing).
- restores aesthetic appearances and marginal integrity when used to repair the interface between restorations and marginal tooth structure.

packaging

4 ml Bottle





α - Bond

The 5th generation bonding system for use in total etching technique

Contains : Different methacrylate monomer & lightcured activators

INDICATIONS FOR USE :

➤ Reliable bonding of LC composite restorations to dental tissues.

packaging

Dropper bottle (5ml)

α – Bond 007

The 7th generation bonding system for use in total self etching technique

Contains : Different methacrylate monomer & lightcured activators

INDICATIONS FOR USE :

➤ Reliable bonding of LC composite restorations to dental tissues.



packaging

Dropper bottle (4ml)

Maxxwhite Office

Professional teeth whitening system

Contains : Hydrogen Peroxide 35 % & comes with desensitizing material

packaging

2 x 0.3 gm power pot 2 x 2.25 liquid
1 x 1.5 gm gingival barrier syringe
2 applicator tips, 2 applicator brush

INDICATIONS FOR USE :

- This is a hydrogen peroxide based in office tooth brightening system requiring minimal chair time.
- It is a neutral pH gel and contains desensitizers to maximize patient comfort.

Maxxwhite Office Plus

Professional teeth whitening system

Contains : 37.5 % Hydrogen Peroxide
Auto mix syringe system contain with
Desensitizer

INDICATIONS FOR USE :

- This is a hydrogen peroxide based in office tooth brightening system requiring minimal chair time.
- It is a neutral pH gel and contains desensitizers to maximize patient comfort.

packaging

6gm Automix Syringe.
Gingival Barrier 1.5gm & Tips



MTA

White | 화이트 | vit | branco | blanc

Reparative cement
복기 시멘트
reparativa cement
cimento reparadora
ciment réparatrice



MTA

White reparative cement

Contains : Calcium Silicates, Calcium aluminates, Calcium aluminoferrite , Calcium Sulfates, radiopacifier

INDICATIONS FOR USE :

➤ External or Internal perforating root resorption treatment preservation of pulp vitality.

packaging

1 gm Powder bottle
5 ml liquid bottle

MTA-JET

Premixed Injectable MTA Paste

Contains : Calcium Silicates, Calcium aluminates, Calcium aluminoferrite, Calcium Sulfates, radiopacifier

INDICATIONS FOR USE :

- MTA Jet is developed for root canal filling and root perforation repair, effective preventing periapical tissue irritation as well as secondary infection.
- Permanent obturation of the root canal system.
- Root Perforation repair.
- Pulp revascularization



packaging

3 x 1gm Syringe



packaging

15 gm powder Jar

10 ml Liquid dropper bottle

AFFIX CEM ONE

Radiopaque glass ionomer luting cement

Contains : Base on fluoro aluminosilicate and polyacrylic acid, tartaric acid

INDICATIONS FOR USE :

- Fixing of crowns , dental bridge work , inlays and onlays Formation of tooth stump.

AFFIX CEM TWO

Radiopaque glass ionomer restorative cement

Contains : Base on fluoro aluminosilicate and polyacrylic acid, tartaric acid

INDICATIONS FOR USE :

- Filling cavities I,II,III and IV glass.
- Art technique filling.
- Primary teeth filling.
- Sealing of fissure small defects , caused by non curious lesion of teeth.
- In aim to from layer base for composite restoration temporary filings.



packaging

15 gm powder Jar

10 ml Liquid dropper bottle



packaging

1 x 2gm porcelain etch
1 x 5 ml porcelain silane
5 ml bond ,5 applicator tips 5
applicator brush.

Porcelain Prep Kit

Contains : 9.6 % hydrofluoric acid
porcelain silane , Bonding Agent

INDICATIONS FOR USE :

➤ Porcelain prep kit is financially advantageous and less invasive to repair , rather than replace a chipped porcelain restoration.



packaging

3 gm Syringe

Natural - Care


Caries removing gel

Composition :

➤ Papain (30K MW), Distiller Water, Methylene blue.

INDICATIONS FOR USE :

➤ Softens demineralized carious dentine prior to mechanical removal during chemo- mechanical minimal invasive treatment of patients diagnosed with caries.



Sparkle

Prophylyxis paste

Contains : Acidic components and abrasive materials

INDICATIONS FOR USE :

- Softening and removal of tartar without damaging of enamel.
- Different flavours available.

packaging

50 gm Jar

Ora Gel

20 % Benzocaine gel

INDICATIONS FOR USE :

➤ Benzocaine Ora Gel: The toothache gel contains 20% benzocaine ora gel which is a highly effective active ingredient that gives precise and instant relief for oral discomfort and toothache pain caused by dentures or sores.

➤ Fast Acting Relief: Using Anadent toothache gel, you'll experience the fast acting relief of toothache pains with the strength of benzocaine. It can also help you to dull the pain long enough for you to fall asleep



packaging

40 gm Jar



packaging

4.5 gm Jar

Lido - Pulp

**Pulp devitalizer paste
(Arsenic free)**

Contains : Paraformaldehyde,
Lidocaine hydrochloride

INDICATIONS FOR USE :

➤ Pulp devitalization before removal or amputation.

Alveolar

(Dry socket paste)

Contains : Dry socket paste 4% Guaicol,
4% eugenol in a Petrolatum base

INDICATIONS FOR USE :

➤ Indicated for the relife of symptoms associated with post extraction alveolar osteitis.

➤ Use Instant pain relife.



packaging

12 gm Jar

ND - Tempfill

(Temporary filling material)

Contains : Based on Zinc oxide , Zinc Sulphate

INDICATIONS FOR USE :

- Temporary tooth filling in caries treatment.
- Set in contact with saliva.



packaging

40 gm Jar

Apex - Seal

Root canal sealing & filing material resin base

Contains : Formin-Titanium Dioxide - 99.9 % pure silver powder – Radiopaque material & excipients – resin liquid

INDICATIONS FOR USE :

- Permanent obturation of root canals of teeth of the secondary dentition with or without the aid of obturation points.



packaging

10 gm powder
8 ml Liquid



packaging

Dropper bottle (5ml)

N –SDF

Silver Diamine Fluoride (38%) Solution

Benefits

- Water : provides a liquid base for the mixture.
- Silver : helps kill bacteria.
- Fluoride : helps your teeth rebuild the materials they ' made of (known as remineralization)
- Ammonia : helps the solution remain concentrated so that it's maximally effective against active resonance.



packaging

Mix Flavour
Qty . 10 X 0.4 ml + Tips

F – Varnish

**5 % Sodium fluoride varnish
with calcium phosphate**

INDICATIONS FOR USE :

- 5% Sodium Fluoride Varnish contains 22,600 ppm fluoride.
- It has a strong desensitizing action when applied to dental surfaces, treating hypersensitivity quickly and easily.
- This product sets rapidly on contact with saliva, resulting in patient comfort and acceptance.
- Varnish will leave a thin film on the teeth after application.



packaging

200ml bottle

Antifog - Solution

Contains : Sodium Sulphosuccinate

INDICATIONS FOR USE :

- Take out require quantity in beaker.
- Deep require item in solution for 10 seconds.
- Take out item form beaker just wipe it with clean cloth.
- The antifog effect will stay for 24 hours .

Roto - Spray

Dental Turbine Spray

INDICATIONS FOR USE :

- Shake well before use.
- Maintenance spray for roto spray turbines , air motors , straight and contra-angle hand pieces..
- Spray for 1 second every mid-day and evening.



packaging

300 ml



packaging

Dropper bottle (5ml)

Caris Detector

Contains : Based on Glycol

INDICATIONS FOR USE :

- To indicate carries location.
- Fast and easy to use.
- Stains only non-reminerizable carious dentin.
- Helps clinicians to preserve healthy tissue.

Endo-Sol

(Gutta Percha solvent)

Contains : xylene base

INDICATIONS FOR USE :

- For used in dissolves gutta percha & zinc-oxide during root canal re-treatment.
- To remove all types of root canal filling materials including Zinc-Oxide-Eugenol cement
- Fitting of the primary gutta percha point in root-canal filling techniques using the lateral condensation method.



packaging

10 ml glass bottle



packaging

Dropper bottle (10ml)

Decens - PF

Contains : Fluoride & Nitrate

INDICATIONS FOR USE :

- For use in post operative and pre operative sensitivity.
- Usage as desensitise exposed root surface in both direct and indirect restoration.

Formo Kresol

Contains : cresol , formalin

INDICATIONS FOR USE :

- Formo kresol is a non specific bacterial medicament most effective against aerobic and anaerobic found in root canals.



packaging

Dropper bottle (15 gm)



packaging

Dropper bottle (15 ml)

Unilik

Universal Liquid For All Gic

INDICATIONS FOR USE :

- Used for all type of glass ionomer cement
- Excellent wetting properties.
- Glass Ionomer formulation.
- Long term clinical data.
- Good colour match to tooth structure.

Prep - Aqua

17 % EDTA Solution

INDICATIONS FOR USE :

- A Chelating Agent for the reaming of root canals.
- Decalcifies canal walls.
- Eases instrumentation of root canals.
- Buffered to a neutral pH.



packaging

Dropper bottle (15 ml)



packaging

10 ml glass bottle

Copal Varnish

Contains : Copal gum in easily volatile solvent

INDICATIONS FOR USE :

- Cavity varnish issued to seal to dentine tubeles to protect acids, saliva and debirs from the pulp.
- It is use under amalgam restrorations to prevent microleakge and under zinc phosphae cement to prevent penetrations of acid the pulp.
- If cavity liners of medicated based care used , varnish is placed after or on top these materials.

Eugenol

INDICATIONS FOR USE :

- Used as an obtundent for the relief of minor dental pain.
- Reacts with zinc oxide to form a sedative dental dressing.



packaging

15 ml glass bottle

ThixoCal

Contains :Calcium Hydroxide



packaging

2 x 3 gm Syringe & dispensing tips

INDICATIONS FOR USE :

- Direct and Indirect pulp capping.
- Apexification and apxogenesis.

Thixo - Prep

Contains :EDTA Salts (17%) and Carbamide Peroxide



packaging

2 x 3 gm Syringe & dispensing tips

INDICATIONS FOR USE :

- Detection and formation of root canal orifice.
- Lubricating dental instruments at chemical and mechanical enlargement of root canals.

Styptic - Gel

(Hemostatic Gel)

Contains : Based on aluminium chloride (25%)



packaging

2 x 3 gm Syringe & dispensing tips

INDICATIONS FOR USE :

- Gingival capillary bleeding .
- Gingival retraction for impression making and cervical caries treatment.
- Treatment of root canals at apical bleeding.

Iodopex

(Calcium Hydroxide root canal filling material)

Contains : Calcium Hydroxide With iodoform (40%)



packaging

1 x 2.2 gm Syringe & dispensing tips

INDICATIONS FOR USE :

- Chronic and all destructive forms of periodontitis.
- Root apexification in case of periodontitis.
- Infection of root perforations and fracture.

Bright - Smile

(Home bleaching gel)

Contains : Carbamide Peroxide (22%)



packaging

2 x 3 gm Syringe & dispensing tips


INDICATIONS FOR USE :

- For tooth Whitening.
- Easy to use.
- With pleasant flavour.

ND-Shine

(Composite polishing paste)

Contains : Aluminium Oxide base



packaging

2 x 2 gm Syringe & dispensing tips

INDICATIONS FOR USE :

- ND Shine polishing paste that creates a beautiful final polish on composites, porcelain, metal and natural tooth structure.

Thixo Etch

(Gel for etching)

Contains : Orthophosphoric (37%)



packaging

2 x 3 gm Syringe & dispensing tips

INDICATIONS FOR USE :

➤ Enamel and dentine etching according to acid etching technique.

Thixo Etch (Jumbo)

(Gel for etching)

Contains : Orthophosphoric (37%)



packaging

1 x 50 gm Syringe & dispensing tips

INDICATIONS FOR USE :

➤ Enamel and dentine etching according to acid etching technique.

ND-Digital X-Ray Sensor Sleeve

- ◆ Made from Food Grade Plastic
- ◆ Softer make for patient comfort
- ◆ Protects sensor contamination
- ◆ No loss of clarity
- ◆ Fits almost all types of sensor

packaging

500 Pcs

Rvg - Sleeve


INDICATIONS FOR USE :

- The quality of Sleeves is ultra-soft that will not hurt the patient tissues.
- Latex Free.
- Generally fits on any sensor.
- Available in long type that also provides cross contamination from sensor wire.
- Easy to use.
- Made of Food Grade Quality.

Zinc –Oxide

INDICATIONS FOR USE :

- 99.00 percent pure.
- Arsenic free, pure zinc oxide powder for temporary dental dressing.
- Makes smooth paste with eugenol.
- Extrafine, odorless, white powder.
- Long shelf life.
- Provides sufficient compressive strength.



packaging

110 gm Bottle

Transforming Healthcare
by driving

A **vailability**
ccessibility
ffordability
wareness

Marketed By:



MEDIKABAZAAR

Boston Ivy Healthcare Solutions Private Limited

3rd Floor, Second Avenue, Atul Projects,
Subhash Nagar, Andheri East, Mumbai,
Maharashtra 400093.



+91 794805 8625



support@medikabazaar.com



www.medikabazaar.com