



ISO 13485 & CE Certified Company

Dealing in Medical Equipments & Services

400W DIGITAL CAUTERY MACHINE

BIOHEARTZ - 400W-P

BIOHEARTZ - 400W-P



ELECTRO SURGICAL UNIT


Technical Specification

		Max.Power(Watts)	RatedLoad(Ohms)
	Operating Frequency	480 Khz	-
	No. of Program	20	-
	Display	7 Segment LED Digital	-
CUT	Pure Cut	400 W	500
	Blend 1	250 W	500
	Blend 2	200 W	500
	Endocut	99 W	500
COAG	Spray	120 W	500
	Fulguration	150 W	500
	Desiccation / Soft	150 W	500
BIPOLAR	Bipolar Coagulation Micro	80 W	100
	Bipolar Coagulation Macro	80 W	100
	Bipolar Cut	99 W	100

Standard Accessories

High quality and safe accessories

Monopolar Double Pedal Footswitch



Bipolar Single Pedal Footswitch




Patient Plate With Cable Cord Autoclavable



Chuk Handle With Electrode Autoclavable



Bipolar Forcep Autoclavable With Cable



Disposable Hand Switch Pencil Imported

