

TANDA

CLASS B AUTOCLAVE



TANDA



ADVANTAGE

- It adopts European B standard with 3-times pre-vacuum preceding vacuum drying. Residual humidity is less than 0.2%.
- Vacuum measurement reaches-0.8bar. It sterilizes and sanitizes miscellaneous packed or unpacked medical instruments, Whether solid or A-type hollow, whether with multi-apertures or inset pipes. It efficiently tides air from every hollow area of any hollow instrument, such as hand-piece, making sure that water steam reaches every comer.
- Mini printer can be attached to record the process of sterilization.



PROGRAM LIST

Program	Temperture	Pressure	Pre-Vacuum	SterilizationTime	DryingTime
UNPACKED	134°C	210kPa	1	4 min	9 min
PACKED	134°C	210kPa	3	5 min	9 min
PRION	134°C	210kPa	3	18 min	9 min
POROUS	121°C	110kPa	3	20 min	18 min

SPECIFICATION

Model	Voltage/Power	Chamber Size	Outside Size	Packing Size	Net Weight
TANDA 23L	220V/1700W	Φ 245X470 mm	675X450X400 mm	770X525X495 mm	53 kg

