

Specifications: A2000



SHENZHEN COMEN MEDICAL INSTRUMENTS CO.,LTD

Floor 10, Floor 11 and Section C of Floor 12 of Building 1A & Floor 1 to
Floor 5 of Building 2, FIYTA Timepiece Building, Nanhuan Avenue,
Matian Sub-district, Guangming District, Shenzhen, Guangdong,
518106,P.R. China

Tel: +86-755-26408879

Fax: +86-755-26431232

Email: info@szcomen.com

Web: www.comen.com

Vaporizer

A2000



Physical Characteristics

Dimensions	210mm*120mm*237mm
Weight	5kg±1kg
Protection level	IP21
Mounting type	Selectatec compatible Dräger plug-in cagemount (ISO 23mm tapers)

Filling Volume

Dry wick	360ml (±60ml)
Moist wick	300ml (±60ml)
Between the MAX and MIN scale marks	260ml (±60ml)

Filler Type

Sevoflurane	Pour Fil
	Quik Fil
	Easy Fil
	Key Fil
Isoflurane	Pour Fil
	Easy Fil
	Key Fil

Performance Specification

Range	Sevoflurane: 0%~8%
	Isoflurane: 0%~6%
Accuracy	≤5% (0.2~15 l/min): ±0.2% or ±20%rel
	≥5% (0.2~10L/min): ±0.2% or ±20%rel
Compensation methods	Pressure compensation: Spiral pipeline
	Flow compensation: Spiral pipeline
	Temperature compensation: Metal valve
Filling speed	≥2ml/s

Operating Environment

Temperature	10~40°C
Humidity	15%~95%
Atmospheric pressure	70.0~106.0kPa
Maximum inclination angle	Stand alone, not fixed: ≤10°
	During use: ≤30°
	In transport: not limited

*Specifications subject to changes without prior notice.

*All rights reserved by Comen.