

Ozocheck Uscan 120

World Class Urine Analyzer



- Compatible with Ozocheck 10P Urine Strips
- In-built thermal printer
- PC/LAN compatible
- Stable in Indian environmental conditions
- 1000 tests Memory
- Prevents cross-infection of urine through inbuilt mechanism



**IDEAL FOR
ALL SIZE
OF LABS**



13485:2016



**MDR
2017**



Easy



Convenient



Reliable



Accurate



Quick

Ozocheck Uscan 120

World Class Urine Analyzer



Feature

Type of Analyser	:	Semi-automated system
Principle	:	Super High luminosity cold light source reflection (Photoelectric Colorimetry)
Compatibility	:	Compatible with Ozocheck Urine Strips 10P
Available Parameters	:	LEU, NIT, URO, PRO, PH, BLD, SG, KET, BIL, GLU
Testing Speed	:	Single Test: 60 strips/hr, Continuously >120 strips/hr
Wavelengths	:	Uses Triple monochromatic wavelengths (558nm, 635nm & 730nm)
Memory	:	1000 tests
Nil Trace Mechanism	:	Absorbs excess urine automatically, thereby preventing cross-inflection
Printer	:	In-built Thermal Printer, with outer stylus printer interface
Connectivity	:	PC connected with RS232
Data Communication	:	Data communication Interface Baud rate of 9600bps
Operating Conditions	:	0-40°C, RH <85%
Language Interface	:	English
Power Requirements	:	AC220V, 50-50Hz, ≤60W
Weight	:	4kg
Dimensions	:	355mm*300mm*145mm



Distributed By:



DIAGNOSOL HEALTHTECH

T-134, Ganesh Nagar, Niwaru Road, Jhotwara, Jaipur- 302012(Raj.)



+91-8432252182



Info.diagnosol@gmail.com



www.diagnosolhealthtech.com