



DEXUP LITHOTOMY STIRRUPS



PRODUCT FEATURES



PADDED FOR COMFORT

Padded leg supports for comfort.



SECURE LEG STRAPS

Straps to secure legs in position.



STABLE DURING SURGERY

Provides stability during surgical procedures.



PRECISE LEG POSITIONING

Allows for precise positioning of legs.



EASY SURGICAL ACCESS

Facilitates access to the surgical site.



SAFE & COMFORTABLE USE

Ensures patient safety and comfort.



ADJUSTABLE FOR ALL PATIENT SIZES

Adjustable design for different patient sizes.








ADVANTAGES

- 1 PROPER LEG POSITIONING**
Facilitates proper positioning of the patient's legs during surgical procedures.
- 2 STABLE & SECURE SUPPORT**
Provides stability and support, minimizing the risk of unintended movement.
- 3 IMPROVED SURGICAL ACCESS**
Allow surgeons to better access to the surgical site, enhancing precision.
- 4 ENHANCED PATIENT COMFORT**
Help keep patient comfort and reduce discomfort during prolonged procedures.
- 5 ENSURES PATIENT SAFETY**
Ensure patient safety by preventing falls or slips during surgery.
- 6 BOOSTS OPERATING ROOM EFFICIENCY**
Enhances overall efficiency in the operating room by aiding in keeping the desired surgical posture.

AREA OF APPLICATIONS

- 1 GYNECOLOGY PROCEDURES**
- 2 UROLOGY PROCEDURES**
- 3 PELVIC SURGERIES**
- 4 LOWER ABDOMINAL SURGERIES**
- 5 PROCEDURES REQUIRING LITHOTOMY POSITION**

CONTRAINDICATIONS

-  Severe peripheral vascular disease.
-  Hip dislocation or instability.
-  Recent lower extremity fractures.
-  Compromised skin integrity or wounds on the legs.
-  Neurological defects affect leg positioning.
-  Uncontrolled hypertension or cardiovascular instability.
-  Known allergies or sensitivity to materials used in the stirrups.

DEVICE SPECIFICATIONS

Cantilever length	1000mm
Applicable patient weight	≤ 180KG
Cantilever up and down adjustment range	+90 °-45 ° (i.e., adjustment angle ≥ 125 °)
Horizontal adjustment range	0 ° to 25 ° and from 0 ° to 30 ° (i.e., the adjustment angle is ≥ 45 °)
The pneumatic power support rod	Original German equipment
Fore and aft adjustment range	330mm

