

ALL SILICONE FOLEY BALLOON CATHETER

**MRI
SAFE**

FIRST IN INDIA

L MED™


MADE FOR COMFORT

L MED
MADE FOR COMFORT


All Silicone Foley Catheter





STERILE
For Single Use Only




LEELAVATHY MEDICAL DEVICES COMPANY


Type	Size(Fr.)	Nominal (Outer Diameter*)(mm)	Balloon Capacity	Catalogue ID
(I) TWO WAY 	6 Fr	2.00	1.5ml, 3ml	LMED001
	8 Fr	2.66	1.5ml, 3ml	LMED002
	10 Fr	3.33	3ml, 5ml	LMED003
	12 Fr	4.0	5ml, 15ml, 5-15ml	LMED004
	14 Fr	4.66	5ml, 15ml, 5-15ml, 20ml	LMED005
	16 Fr	5.33	5ml, 15ml, 5-15ml, 20ml	LMED006
	18 Fr	6.0	5ml, 15ml, 5-15ml, 20ml	LMED007
	20 Fr	6.66	5ml, 15ml, 5-15ml, 20ml	LMED008
	22 Fr	7.33	5ml, 15ml, 5-15ml, 20ml	LMED009
	24 Fr	8.0	5ml, 15ml, 5-15ml, 20ml	LMED010


Type	Size(Fr.)	Nominal (Outer Diameter*)(mm)	Balloon Capacity	Catalogue ID
(II) THREE WAY 	12 Fr	4.0	5ml, 15ml, 5-15ml	LMED011
	14 Fr	4.66	5ml, 15ml, 5-15ml, 20ml	LMED012
	16 Fr	5.33	5ml, 15ml, 5-15ml, 20ml	LMED013
	18 Fr	6.0	5ml, 15ml, 5-15ml, 20ml	LMED014
	20 Fr	6.66	5ml, 15ml, 5-15ml, 20ml	LMED015
	22 Fr	7.33	5ml, 15ml, 5-15ml, 20ml	LMED016
	24 Fr	8.0	5ml, 15ml, 5-15ml, 20ml	LMED017

Type	Size(Fr.)	Nominal (Outer Diameter*)(mm)	Balloon Capacity	Catalogue ID
(III) TWO WAY WITH CAP 	6 Fr	2.00	1.5ml, 3ml	LMED018
	8 Fr	2.66	1.5ml, 3ml	LMED019
	10 Fr	3.33	3ml, 5ml	LMED020
	12 Fr	4.0	5ml, 15ml, 5-15ml	LMED021
	14 Fr	4.66	5ml, 15ml, 5-15ml, 20ml	LMED022
	16 Fr	5.33	5ml, 15ml, 5-15ml, 20ml	LMED023
	18 Fr	6.0	5ml, 15ml, 5-15ml, 20ml	LMED024
	20 Fr	6.66	5ml, 15ml, 5-15ml, 20ml	LMED025
	22 Fr	7.33	5ml, 15ml, 5-15ml, 20ml	LMED026
	24 Fr	8.0	5ml, 15ml, 5-15ml, 20ml	LMED027


Type	Size(Fr.)	Nominal (Outer Diameter*)(mm)	Balloon Capacity	Catalogue ID
(IV) THREE WAY WITH CAP 	12 Fr	4.0	5ml, 15ml, 5-15ml	LMED028
	14 Fr	4.66	5ml, 15ml, 5-15ml, 20ml	LMED029
	16 Fr	5.33	5ml, 15ml, 5-15ml, 20ml	LMED030
	18 Fr	6.0	5ml, 15ml, 5-15ml, 20ml	LMED031
	20 Fr	6.66	5ml, 15ml, 5-15ml, 20ml	LMED032
	22 Fr	7.33	5ml, 15ml, 5-15ml, 20ml	LMED033
	24 Fr	8.0	5ml, 15ml, 5-15ml, 20ml	LMED034


SUPRA PUBIC CATHETER SPECIFICATION

Type	Size(Fr.)	Nominal (Outer Diameter*)(mm)	Balloon Capacity	Catalogue ID
(I) TWO WAY 	10 Fr	3.33	3ml, 5ml	LMED035
	12 Fr	4.0	5ml, 15ml, 5-15ml	LMED036
	14 Fr	4.66	5ml, 15ml, 5-15ml, 20ml	LMED037
	16 Fr	5.33	5ml, 15ml, 5-15ml, 20ml	LMED038
	18 Fr	6.0	5ml, 15ml, 5-15ml, 20ml	LMED039
	20 Fr	6.66	5ml, 15ml, 5-15ml, 20ml	LMED040
	22 Fr	7.33	5ml, 15ml, 5-15ml, 20ml	LMED041
	24 Fr	8.0	5ml, 15ml, 5-15ml, 20ml	LMED042

Type	Size(Fr.)	Nominal (Outer Diameter*)(mm)	Balloon Capacity	Catalogue ID
(II) THREE WAY 	12 Fr	4.0	5ml, 15ml, 5-15ml	LMED043
	14 Fr	4.66	5ml, 15ml, 5-15ml, 20ml	LMED044
	16 Fr	5.33	5ml, 15ml, 5-15ml, 20ml	LMED045
	18 Fr	6.0	5ml, 15ml, 5-15ml, 20ml	LMED046
	20 Fr	6.66	5ml, 15ml, 5-15ml, 20ml	LMED047
	22 Fr	7.33	5ml, 15ml, 5-15ml, 20ml	LMED048
	24 Fr	8.0	5ml, 15ml, 5-15ml, 20ml	LMED049


TIP HOLE FOLEY CATHETER SPECIFICATION


Type	Size(Fr.)	Nominal (Outer Diameter*)(mm)	Balloon Capacity	Catalogue ID
(I) TWO WAY 	12 Fr	4.0	5ml, 15ml, 5-15ml	LMED050
	14 Fr	4.66	5ml, 15ml, 5-15ml, 20ml	LMED051
	16 Fr	5.33	5ml, 15ml, 5-15ml, 20ml	LMED052
	18 Fr	6.0	5ml, 15ml, 5-15ml, 20ml	LMED053
	20 Fr	6.66	5ml, 15ml, 5-15ml, 20ml	LMED054
	22 Fr	7.33	5ml, 15ml, 5-15ml, 20ml	LMED055
	24 Fr	8.0	5ml, 15ml, 5-15ml, 20ml	LMED056

Type	Size(Fr.)	Nominal (Outer Diameter*)(mm)	Balloon Capacity	Catalogue ID
(II) THREE WAY 	12 Fr	4.0	5ml, 15ml, 5-15ml	LMED057
	14 Fr	4.66	5ml, 15ml, 5-15ml, 20ml	LMED058
	16 Fr	5.33	5ml, 15ml, 5-15ml, 20ml	LMED059
	18 Fr	6.0	5ml, 15ml, 5-15ml, 20ml	LMED060
	20 Fr	6.66	5ml, 15ml, 5-15ml, 20ml	LMED061
	22 Fr	7.33	5ml, 15ml, 5-15ml, 20ml	LMED062
	24 Fr	8.0	5ml, 15ml, 5-15ml, 20ml	LMED063

Note : Available In Standard & Custom Length


TIEMANN OR COUDE TIP CATHETER

Type	Size(Fr.)	Nominal (Outer) Diameter*(mm)	Balloon Capacity	Catalogue ID
(I) TWO WAY 	12 Fr	4.0	5ml, 15ml, 5-15ml	LMED064
	14 Fr	4.66	5ml, 15ml, 5-15ml, 20ml	LMED065
	16 Fr	5.33	5ml, 15ml, 5-15ml, 20ml	LMED066
	18 Fr	6.0	5ml, 15ml, 5-15ml, 20ml	LMED067
	20 Fr	6.66	5ml, 15ml, 5-15ml, 20ml	LMED068
	22 Fr	7.33	5ml, 15ml, 5-15ml, 20ml	LMED069
	24 Fr	8.0	5ml, 15ml, 5-15ml, 20ml	LMED070

Type	Size(Fr.)	Nominal (Outer) Diameter*(mm)	Balloon Capacity	Catalogue ID
(II) THREE WAY 	16 Fr	5.33	5ml, 15ml, 5-15ml, 20ml	LMED071
	18 Fr	6.0	5ml, 15ml, 5-15ml, 20ml	LMED072
	20 Fr	6.66	5ml, 15ml, 5-15ml, 20ml	LMED073
	22 Fr	7.33	5ml, 15ml, 5-15ml, 20ml	LMED074
	24 Fr	8.0	5ml, 15ml, 5-15ml, 20ml	LMED075

SILICON TUBES

Specification : ID 0.2mm to ID 12mm Length : Custom Made



SPECIAL FEATURES

- Radio opaque tip for clear identification of the catheter position
- Additional silicone coating to ensure extra smooth penetration
- Uniform roundness throughout the length of the catheter tube.
- Smooth tip facilitates for painless insertion.
- Balloon integration with the catheter is done with a special process to eliminate rough edges.
- Comparitively bigger drainage lumen for better and increased Drainage flow than in standard latex catheter.
- Uniform and larger inflation lumen for easy inflation and deflation of balloon.
- 100% biocompatible silicone is used and therefore biological activity is nil
- Clear silicone material facilitates easy visibility of drainage.

REFERENCE DOCUMENTS

This Product conforms to - BS EN ISO 20696 : 2018

RECOMMENDED STORAGE CONDITION


Store at room temperature (30-35^oc / 86-95^of) and avoid exposure To elevated temperatures. Do not expose product to radiation, Sunlight and uv light, as it may change the properties of the product

PACKING CONFIGURATION


10 Catheter in 1 Box (Unit Carton Box), 10 Unit Carton in 1 master Carton (Shipping carton Box)

IN-HOUSE TESTING


EYE OCCULTATION




TIP INTEGRITY




WATER RECOVERY TEST



FLOW RATE TEST



TUBE INTEGRITY





**Class 10000
Clean Room
Facility**



LEELAVATHY MEDICAL DEVICES COMPANY

No.1/122, Paraniputhur-Kovur Service Road, Periyapanicherry,
Chennai-600128, INDIA

Email Id : sales@lmed.co.in | Web : www.lmed.co.in

Phone: +91 8939480062, 9884399735



ISO 13485 Certified
WHO GMP Certified

