



## **CoaguChek® XS Pro System**

*Offering barcode scanning technology and options for connectivity*



## The CoaguChek XS Pro System

*The CoaguChek XS Pro system ensures accurate PT/INR results and offers data management capabilities for healthcare professionals who manage anticoagulation patients at the point of care.*

### Accurate and flexible technology

- Accurate, precise results in about 1 minute
- Flexible QC options for any testing environment: onboard controls and optional liquid quality controls
- Only handheld, PT/INR monitoring system that neutralizes therapeutic levels of heparin (up to .8 U/mL) and LMWH (up to 2 IU/mL)
- INR corrected for wide range of hematocrit<sup>1</sup>

### Streamlined data management

- Stores 2000 patient test results and 500 optional liquid quality control results
- Operator and LQC Lockout capabilities (\*Operator lockout only available when using with data management system)
- Stores up to 60 code keys
- Ability to denote patient ID in several fields
- Ability to add comments after each test result

### Easy and convenient to use


- Small sample size (8 µL) obtained by finger stick
- Runs on 4 AA batteries or optional rechargeable battery pack
- Icon-driven, color, touch screen interface

### Integrated barcode scanner

- Provides quick, simple way to capture patient IDs and operator IDs
- Virtually eliminates manual entry errors

### The One and Only CoaguChek brand: the market leader in point-of-care anticoagulation monitoring

Since 1994, Roche Diagnostics has helped healthcare professionals provide quality care to their patients with a PT/INR device they trust. Find out why more CoaguChek test strips are sold for point-of-care anticoagulation testing than all other brands combined.<sup>2</sup>



CoaguChek XS Pro PT Test Strip package insert, 05967716001(02), 2010.