

Runyes[®]
FOCUS ON DENTAL

Ray 98 (P)

Camera Type



Gun Type



Portable High-Frequency X-Ray Unit

Compact & Wireless DC X-Ray

Camera Type

Technical Specifications	
Focal Spot	0.4 mm
Tube Current	2mA + 20%
Voltage	70 kV + 10%
Frequency	50 -60Hz
Total Filtration	1.5 AL /70kv
Exposure Time Range	0.04 Sec – 2.0Sec
Extra Cone Coating	0.5 mm



Available for hand-held shooting separately

- ▲ Useful for multiple chair Clinical Set-up
- ▲ Very Light in weight device, easy for single handed operation with the weight of 1.7 Kg
- ▲ Superior Image Quality with less radiation
- ▲ Imported circuits from South Korea, ensures high-quality image even during low battery power
- ▲ 150 shots of X-Ray can be taken once the unit is fully charged
- ▲ Various exposure modes for different parameters
- ▲ Wide digital display for easy operation
- ▲ Powered by 12v /2500mAH lithium ion battery backup
- ▲ Complete body lead coating of 0.8mm
- ▲ Total emissivity-1 mtr position<0.25mGrey/h(70kv /2mA/1s loading)



Easy Operation



More applicable



Clear image



Portable use



Concise & clear digital display

Handheld

Portable X-Ray

IMAGING IN MY CONTROL

MODERN MEDICAL IMAGING TECHNOLOGY
 More accurate and more flexible.

The Gun-Type Portable X-Ray machine is a cutting-edge choice for dental diagnosis, featuring ergonomic design, long battery life, and intelligent angle adjustment for precise imaging.

Gun Type

Smart angle adjustment
 for clearer images



Large screen clearly display

Button adjustment to improve
 diagnostic efficiency

Long handle, comfortable to
 hold, easier to operate



Atmospheric & concise design

Recessed power button. Avoid accidental touch



Smart camera
 angle assistance



Ergonomic design



Long battery life

Technical Specifications

Tube Voltage	70kV±10%
Tube Current	2mA + 20%
Input Voltage	100V-240V
Frequency	50 -60Hz
Anode Angle	12.5 °
Focus	0.4mm
Exposure Time Range	0.04s -2.0s
Leak Emissivity	1 meter 0.25mGy/h
Machine Weight	1.7kg
Relative Humidity	30%RH-75%RH

Transforming Healthcare by driving

A
availability
accessibility
affordability
awareness

Marketed By:



MEDIKABAZAAR

Boston Ivy Healthcare Solutions Private Limited

3rd Floor, Second Avenue, Atul Projects,
Subhash Nagar, Andheri East, Mumbai,
Maharashtra 400093.



+91 794805 8625



support@medikabazaar.com



www.medikabazaar.com