

Healix Conmem PRO GOLD

Variable Taper Rotary Files



- ▶ A-R Wire Tech for Shaping Files
- ▶ A-H Wire Tech for Finishing Files
- ▶ 3 Shaping Files, 3 Finishing Files
- ▶ Specialty for Variable Taper 6 Files Design
- ▶ Suitable for Curved Root Canal Preparation
- ▶ High elasticity of NITI Material

6 PCS/BOX

# ISO	% Taper	File Colour	Length Available	
SX	8%	Gold	21mm, 25mm	
S1	2%	Purple	21mm, 25mm	
S2	4%	White	21mm, 25mm	
F1	6%	Yellow	21mm, 25mm	
F2	7%	Red	21mm, 25mm	
F3	8%	Blue	21mm, 25mm	

6 PCS/BOX

Assorted		Sx, S1, S2, F1, F2, F3	21mm, 25mm	
----------	--	------------------------	------------	--



Healix Conmem M3 Gold

Constant Taper Rotary Files



Conmem Gold Wire technology

Extremely flexible files

Files of choice for very narrow and calcified canals

Files of choice for severely curved canals

Minimal invasive files

Improved resistance to cyclic fatigue

TOP PRECISION AUTO MANUFACTURING

Advanced control memory alloy tech
Conmem Gold wire Improved resistant
Much better cyclic fatigue resistance
Reserve more dentine

INACTIVE TIP

Advanced tip process to avoid forming steps
Suitable for curved canals

ENDO STOP

Mark the orient of the root
canal and the working length

6 PCS/BOX

# ISO	% Taper	File Colour	Length Available	
#17	4%	White	21mm, 25mm	○
#20	4%	Yellow	21mm, 25mm	●
#25	4%	Red	21mm, 25mm	●
#30	4%	Blue	21mm, 25mm	●
#17	6%	Double White	21mm, 25mm	○ ○
#20	6%	Double Yellow	21mm, 25mm	● ●
#25	6%	Double Red	21mm, 25mm	● ●
#30	6%	Double Blue	21mm, 25mm	● ●

6 PCS/BOX

Assorted	4%,6%	Orifice Opener, #17-4%, #20,4% #20-6%, #25-4% #25-6%,	21mm, 25mm	● ○ ● ● ● ●
----------	-------	--	------------	----------------------------