

ReeLife™

Video Laryngoscope



CATALOGUE



GeM
Government
e Marketplace

Website : www.reliefsurgico.in

ReeLife™

Video Laryngoscope



Key Features & Advantages

- **Craftsmanship**
 - Ergonomic design, easy to hold.
 - Durable and high quality materials are adopted to prevent breakage
- **Clear Horizon**
 - 2 MP camera & dual antilog lens
 - 2.8" color LCD, anti-reflection display
 - Adjust screen angle, up / down $\geq 120^\circ$ more convenient to observe
- **One Button Operation**
 - Simple to operate
 - Continuous oper time: > 3 hours
 - The disposable blade can be replaced quickly and conveniently

