



India's Most Trusted Suction Brand

ISO
9001:2015

ISO
13485:2016



FLOVAC 50

Main Applications

- Emergency Department
- Intensive Care Unit
- OT / Surgery
- Obstetrics
- Gastroenterology

ELECTRONIC SENSOR



Convenient emptying
handles on Lid & Jar
- Controls Infection

JAR



TECHNICAL FEATURES

Model	HI VAC PLUSS VARIABLE
Housing	Plastic moulded on two wheels with foot switch & changeover lever
Max. Vacuum / LPM	700 mm Hg, 60 Ltrs / min
Pump Type	piston pump
Jars	2 x 4 Ltrs. Polysulfone Jars
Filter	Bacterial Filter Autoclavable / Reusable
Tubing	10 mm ID x 2 mtr. (PVC)
Vacuum Gauge	6.25 cm dia 0-760 mm Hg
Power Supply	220 V AC, 50 / 60 Hz, 180 W
Noise Level	55 dB A $\pm$ 3
Dimensions	46 x 34 x 86 cms
Net Weight	22 Kg
Overflow Safety	Mechanical & Electronic Sensor
Optional	110 V
Certifications	CE, IEC 60601-1-2012