

NEOALGIN

Alginate Impression Material



Colour Change Pattern: Non-Chromatic
(Sky Blue, No colour change)

Powder:Water Ratios: 18gm:36ml
(Warm water and hard water accelerates the setting time, Cold water retards the setting time)

Impression Count: 25 Impressions
(Based on 18gm needed for one Impression of tray 2/3 or medium tray size)

Features

- Sky blue colour
- Spearmint flavor for high acceptance
- Mixing time: 60 Sec
- Time in mouth: 30 Sec
- Total working time: 60 Sec